



## Dean Madden

National Centre for Biotechnology Education, University of Reading  
Science and Technology Centre, Earley Gate, Reading RG6 6BZ UK | E: D.R.Madden@reading.ac.uk

# Better milk for cats

## Immobilised lactase used to make lactose-reduced milk

### Aims

This simple practical investigation introduces students to the principles of digestion and enzyme immobilisation. It can be used as the starting point for numerous more advanced activities such as the regulation of lactase production in *Escherichia coli* (the *lac* operon), the evolution and social significance of lactose tolerance in humans and the use of enzymes in food production.

### Introduction

Lactase ( $\beta$ -galactosidase) catalyses the hydrolysis of lactose to glucose and galactose. Both of these sugars taste sweeter than lactose and are more readily-digestible than them. Despite their traditional liking for milk, cats are unable to digest large amounts of lactose. Milk can be treated with the enzyme to make a lactose-reduced milk suitable for cats or humans who are lactose intolerant.

Although the production of a special 'cat milk' may seem trivial, an estimated 75% of the world's population are intolerant of lactose in adulthood — it is lactose tolerance that is unusual.

Commercially, the milk is treated by injecting enzyme into the carton as UHT milk is packaged, or by using an immobilised enzyme — an enzyme which has been trapped on an inert material so that it can be used repeatedly.

In this activity, students immobilise the lactase in calcium alginate beads held within a small column, over which the milk is passed.



## Equipment and materials

### Needed by each person or group

#### Equipment

- Small piece (about 1 cm<sup>2</sup>) of nylon gauze *e.g.*, net curtain
- 10 mL plastic syringe (without a needle)
- 4 mm diameter aquarium airline or silicone tubing, about 7 cm long, to fit syringe
- Aquarium airline tap or adjustable laboratory tubing clip (Hoffman clip)
- Retort stand, boss and clamp (to support enzyme column)
- Small beakers or disposable plastic cups, 2
- Tea strainer
- Glass stirring rod

#### Materials

- Lactase enzyme, *Novozymes Lactozym*®, 2 mL
- 2% sodium alginate solution, 8 mL
- 1.5% calcium chloride solution, 100 mL
- Milk, 50 mL, ideally at room temperature (*not* UHT milk)
- Semi-quantitative glucose test strips  
*e.g.*, *Roche Diabur-Test 5000* or *Ames Diastix*

#### Note

*All solutions must be made up using distilled or deionised water. Sodium alginate is not readily soluble, and requires both warm water and stirring to dissolve it.*

## Procedure

- 1 Mix the enzyme with the sodium alginate solution, then draw it up into a 10 mL syringe.
- 2 Add the alginate-enzyme mixture a drop at a time from the syringe to the calcium chloride solution. *Do not allow the tip of the syringe to come into contact with the calcium chloride solution, as this will cause the alginate to harden, blocking the outlet.* The beads, which contain the enzyme immobilised in a matrix of calcium alginate, should be allowed to harden for a few minutes.
- 3 Attach a short length of tubing to the tip of a syringe barrel. Place a small disc of nylon gauze inside the barrel, to prevent the beads from blocking the syringe outlet.

Fig. 1




Fig. 2





Fig. 3



- 4 Separate the beads of immobilised enzyme from the liquid with the tea strainer.
- 5 Carefully tip the beads into the syringe barrel.
- 6 Close the tubing on the syringe barrel using a tubing clip.

Fig. 4




Fig. 5





Fig. 6



- 7 Test the milk before treatment using the glucose test strips, to ensure that it does not contain any glucose.
- 8 Pour a small volume of milk over the enzyme beads, then undo the clip and allow the treated milk to run into a small beaker.
- 9 Test the milk leaving the column using the glucose test strips. If necessary, return the treated milk to the column until the desired concentration of glucose is achieved.

Fig. 7




Fig. 8





Fig. 9



## Safety guidelines



### Do not consume the milk

The enzyme suggested for this work is safe to use, provided it is handled appropriately. While *Novozymes Lactozym*<sup>®</sup> is a food-grade products, milk prepared using it should not be consumed. This is because the enzyme has not been handled aseptically, so it (and the product made using it) may have been contaminated.

Readers are advised to refer to any local safety guidelines and to carry out their own risk assessment for any practical work.

### General enzyme safety guidelines

As enzymes are water-soluble, water should always be used for their removal if they are spilt.

### Do not let liquid enzyme preparations dry up

If liquid preparations are allowed to dry up, there is a risk of dust formation. In susceptible people the repeated inhalation of such dust may provoke asthma or a reaction similar to hay fever. Any spillage — on equipment, on the floor or bench — should immediately be rinsed away with water.

### Avoid the formation of aerosols

If enzyme-containing aerosols are formed, there is a risk of inhalation of the enzyme. In susceptible people the repeated inhalation of such aerosols may provoke asthma or hay fever. For this reason enzyme preparations should never be sprayed.

### Avoid direct skin and eye contact

If, by accident, you get liquid enzyme on your skin or in your eyes, the remedy is plenty of tap water. The same applies to clothing. In the event of a spill on clothes, rinse with water then wash as usual. This treatment will generally prove sufficient, but if symptoms develop in the respiratory passages, on the skin or in the eyes, consult a doctor immediately.

## Preparation and timing

This activity takes about 40 minutes. The sodium alginate takes some time to dissolve, so the solution is best prepared before the lesson. The immobilised enzyme may be prepared in advance if desired: the beads should be refrigerated, although they will not keep for more than a few days.

## Troubleshooting

Some UHT milk will test positive for glucose, probably because the heat treatment hydrolyses some of the lactose. UHT milk should therefore be avoided.

## Additional investigations

- 1 The immobilised enzyme column may also be used to treat whey, producing a sweet whey syrup, which is widely used in confectionery (it is usually described on labels as 'hydrolysed whey syrup' or just 'whey syrup').



*Lactose-reduced milk for cats*



*This Italian low-lactose milk is produced using immobilised  $\beta$ -galactosidase to treat semi-skimmed milk*

- 2 Lactase is strongly inhibited by galactose (one of the products of its action on lactose). Hence the flow rate of the substrate over the column is critical to the rate of the enzyme-catalysed reaction: too fast and there isn't time for the reaction to occur; too slow a rate and galactose will accumulate and then inhibit the reaction. Students can therefore investigate the effect of flow rate on the conversion of lactose to glucose and galactose.

## Other sources of information

### Publications

*Immobilised enzymes and cells: a practical approach* by Jonathan Woodward [Ed] (1985) Oxford University Press, Oxford. ISBN: 0 947946 21 7.

*An academic laboratory manual describing methods of immobilising enzymes and cells.*

Richmond, M.L., J.I. Gray and C.M. Stine (1981) Beta-galactosidase: Review of recent research related to technological application, nutritional concerns and immobilization. *Journal of Dairy Science* **64** (9) 1759–1771.

Bayless, T.M., D.M. Paige and G.D. Ferry (1971) Lactose intolerance and milk drinking habits *Gastroenterology* **60** (4) 605–608.

### Web sites

Novozymes A/S, Denmark  
<http://www.novozymes.com>

Wikipedia: Lactose intolerance  
[http://en.wikipedia.org/wiki/Lactose\\_intolerance](http://en.wikipedia.org/wiki/Lactose_intolerance)

## Suppliers

The NCBE supplies *Novozymes* enzyme products to schools and colleges in the UK. Similar arrangements may exist in other countries.

## Storage of materials

The enzyme preparation should be stored, undiluted, at 3–4 °C.

## Acknowledgement



This practical protocol was adapted for the *Volvox* project, which is funded under the Sixth Framework Programme of the European Commission.